

Lot No.: 5H8

## Olerup SSP<sup>®</sup> DRB1\*04 Add-on

Product number: 101.814-12 – including *Taq* polymerase  
101.814-12u – without *Taq* polymerase  
Lot number: 5H8  
Expiry date: 2023-03-01  
Number of tests: 12  
Number of wells per test: 3+1

### CHANGES COMPARED TO THE PREVIOUS HLA-DRB1\*04 ADD-ON LOT (4F1):

The DRB1\*04 Add-on primer set is unchanged compared to the previous *Olerup SSP<sup>®</sup> DRB1\*04 Add-on (Lot No. 4F1)*.

THE NUMBER OF WELLS is unchanged.

### ALLELE COVERAGE:

The DRB1\*04:07 and DRB1\*04:92 alleles give different patterns in the DRB1\*04 Add-on subtyping kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2019-January-23, release 3.35.0.

The primer set also resolves the SBT heterozygous ambiguities:

DRB1\*04:02:01, 04:03:01 = DRB1\*04:37, 04:88

DRB1\*04:02:01, 04:04:01 = DRB1\*04:37, 04:56

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

### RESOLUTION IN HOMO- AND HETEROZYGOTES:

Excellent.

### INFLUENCE ON THE INTERPRETATION OF DRB1\*04 ADD-ON SUBTYPINGS BY NON-DRB1\*04 ALLELES:

None.

Changes in revision R01 compared to R00:

1. The expiration date has been altered due to shelf-life extension.