

Olerup SSP[®] DRB1*04 Add-on

Product number:	101.814-12 – including <i>Taq</i> polymerase 101.814-12u – without <i>Taq</i> polymerase
Lot number:	5H8
Expiry date:	2021-09-01
Number of tests:	12
Number of wells per test:	3+1

The DRB1*04 Add-on primer set is unchanged compared to the previous *Olerup SSP[®] DRB1*04 Add-on (Lot No. 4F1)*.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The DRB1*04:07 and DRB1*04:92 alleles give different patterns in the DRB1*04 Add-on subtyping kit¹; www.ebi.ac.uk/imgt/hla, 2019-January-23, release 3.35.0.

The primer set also resolves the SBT heterozygous ambiguities:

DRB1*04:02:01, 04:03:01 = DRB1*04:37, 04:88

DRB1*04:02:01, 04:04:01 = DRB1*04:37, 04:56

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DRB1*04 ADD-ON SUBTYPINGS BY NON-DRB1*04 ALLELES:

None.