

DRB1\*07  
101.118-24 – licensed for PCR  
101.118-24u – not licensed for PCR  
Lot No.: **X68**

1

## Release Note - *Olerup* SSP™ DRB1\*07 – Lot No. **X68**

Product numbers: 101.118-24 – licensed for PCR  
101.118-24u – not licensed for PCR  
Lot number: X68  
Expiry date: 2009-April-01  
Number of tests: 24  
Number of tubes per test: 8

### CHANGES COMPARED TO THE PREVIOUS DRB1\*07 LOT (V66):

Tube	5'-primer	3'-primer	rationale
3	-	Added	Primer added for the DRB1*0712 allele.
8	-	Added	Primer added for the DRB1*0712 allele.

THE NUMBER OF WELLS is unchanged.

### ALLELE COVERAGE:

DRB1\*0701 to DRB1\*0712, i.e. all the currently recognized DRB1\*07 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2007-April-12, release 2.17.0.

### RESOLUTION IN DRB1\*07 HOMOZYGOTES:

The DRB1\*0701,0701 genotype gives rise to a unique amplification pattern.

### INFLUENCE ON THE INTERPRETATION OF DRB1\*07 SUBTYPINGS BY NON-DRB1\*07 ALLELES:

None.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.