

Olerup SSP[®] DRB1*07

Product numbers: 101.118-24 – including *Taq* polymerase
101.118-24u – without *Taq* polymerase
Lot number: 67F
Expiry date: 2011-April-01
Number of tests: 24
Number of wells per test: 14

CHANGES COMPARED TO THE PREVIOUS DRB1*07 LOT (32E):

Well	5'-primer	3'-primer	rationale
4	-	Modified	Increased specificity of specific primer pair, exchanged positive control primer pair.
13	New	New	Primer pair for the DRB1*0716 allele.
14	New	New	Primer pair for the DRB1*0717 allele.

THE NUMBER OF WELLS has been increased from 12 to 14.

ALLELE COVERAGE:

DRB1*0701 to DRB1*0717, i.e. all the currently recognized DRB1*07 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2009-April-17, release 2.25.2.

RESOLUTION IN DRB1*07 HOMOZYGOTES:

The DRB1*0701,0701 genotype gives rise to a unique amplification pattern.

INFLUENCE ON THE INTERPRETATION OF DRB1*07 SUBTYPINGS BY NON-DRB1*07 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.