

Olerup SSP[®] DRB1*07

Product numbers: 101.118-24 – including *Taq* polymerase
101.118-24u – without *Taq* polymerase
Lot number: 26M
Expiry date: 2012-September-01
Number of tests: 24
Number of wells per test: 16

CHANGES COMPARED TO THE PREVIOUS DRB1*07 LOT (29K):

Well	5'-primer	3'-primer	rationale
10	-	Added	3'-primer added for the DRB1*07:18 allele.
12	Added		5'-primer added for the DRB1*07:21 allele.
13	Added		5'-primer added for the DRB1*07:21 allele.
16	New	New	New primer pairs for the DRB1*07:18 and DRB1*07:20 alleles.

THE NUMBER OF WELLS has been increased from 15 to 16.

ALLELE COVERAGE:

DRB1*07:01 to DRB1*07:21, i.e. all the currently recognized DRB1*07 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2011-April-08, release 3.4.0.

RESOLUTION IN DRB1*07 HOMOZYGOTES:

The DRB1*07:01,07:01 genotype gives rise to a unique amplification pattern.

INFLUENCE ON THE INTERPRETATION OF DRB1*07 SUBTYPINGS BY NON-DRB1*07 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.