

Olerup SSP[®] DRB1*08

Product number: 101.127-12/04 – including *Taq* pol.
101.127-12u/04u – without *Taq* pol.
Lot number: 72E
Expiry date: 2010-June-01
Number of tests: 12 tests – Product No. 101.127-12
4 tests – Product No. 101.127-04
Number of wells per test: 24

CHANGES COMPARED TO THE PREVIOUS DRB1*08 LOT (X72):

Well	5'-primer	3'-primer	rationale
13	-	Added	Primer added for the DRB1*0833 allele.
22	-	Added	Primer added for the DRB1*0834 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DRB1*0801 to DRB1*0834, i.e. all the currently recognized DRB1*08 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-April-08, release 2.21.0.

RESOLUTION IN DRB1*08 HOMOZYGOTES:

The DRB1*0801,0801, DRB1*0802,0802 and DRB1*0804,0804 genotypes give rise to unique amplification patterns.

INFLUENCE ON THE INTERPRETATION OF DRB1*08 SUBTYPINGS BY NON-DRB1*08 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.