

DRB1*09
101.120-06 – licensed for PCR
101.120-06u – not licensed for PCR
Lot No.: Y45

1

Release Note - Olerup SSP™ DRB1*09 – Lot No. Y45

Product number: 101.120-06 – licensed for PCR
101.120-06u – not licensed for PCR
Lot number: Y45
Expiry date: 2009-September-01
Number of tests: 6
Number of tubes per test: 7

Changes compared to the previous DRB1*09 Lot (V79):

Tube	5'-primer	3'-primer	rationale
7	New	New	New primer pair for the DRB1*0906 allele.

THE NUMBER OF WELLS has been increased from 6 to 7.

ALLELE COVERAGE:

DRB1*0901 to DRB1*0906, i.e. all the currently recognized DRB1*09 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2007-October-05, release 2.19.0.

RESOLUTION IN DRB1*09 HOMOZYGOTES:

Very good.

INFLUENCE ON THE INTERPRETATION OF DRB1*09 SUBTYPINGS BY NON-DRB1*09 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.