

## Olerup SSP<sup>®</sup> DRB1\*09

**Product number:** 101.120-06 – including *Taq* polymerase  
101.120-06u – without *Taq* polymerase  
**Lot number:** 18F  
**Expiry date:** 2010-October-01  
**Number of tests:** 6 tests – Product No. 101.120-06  
**Number of wells per test:** 8

### Changes compared to the previous DRB1\*09 Lot (Y45):

Well	5'-primer	3'-primer	rationale
8	New	New	New primer pair for the DRB1*0907 allele.

**THE NUMBER OF WELLS** has been increased from 7 to 8.

### ALLELE COVERAGE:

DRB1\*0901 to DRB1\*0907, i.e. all the currently recognized DRB1\*09 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2008-July-11, release 2.22.0.

### RESOLUTION IN DRB1\*09 HOMOZYGOTES:

Very good.

### INFLUENCE ON THE INTERPRETATION OF DRB1\*09 SUBTYPINGS BY NON-DRB1\*09 ALLELES:

None.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.