

101.120-06 – including *Taq* polymerase101.120-06u – without *Taq* polymeraseLot No.: **60K****Olerup SSP® DRB1*09**

Product number: 101.120-06 – including *Taq* polymerase
101.120-06u – without *Taq* polymerase
Lot number: 60K
Expiry date: 2013-April-01
Number of tests: 6
Number of wells per test: 8

Changes compared to the previous DRB1*09 Lot (07G):

Well	5'-primer	3'-primer	rationale
6	-	Added	New primer added for the DRB1*09:09 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DRB1*09:01 to DRB1*09:09, i.e. all the currently recognized DRB1*09 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-July-16, release 3.1.0.

RESOLUTION IN DRB1*09 HOMOZYGOTES:

Very good.

INFLUENCE ON THE INTERPRETATION OF DRB1*09 SUBTYPINGS BY NON-DRB1*09 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.