

DRA
101.131-24 – licensed for PCR
101.131-24u – not licensed for PCR
Lot No.: **Y51**

1

Release Note - Olerup SSP™ DRA – Lot No. Y51

Product number:	101.131-24 – licensed for PCR 101.131-24u – <u>not</u> licensed for PCR
Lot number:	Y51
Expiry date:	2009-September-01
Number of tests:	24
Number of tubes per test:	2

CHANGES COMPARED TO THE PREVIOUS DRA LOT (R94):

The DRA primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DRA*0101 to DRA*0102, i.e. all the currently recognized DRA alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2007-October-05, release 2.19.0.

RESOLUTION IN DRA HOMOZYGOTES:

Excellent.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.