

Olerup SSP® DRB1*14:01/14:54

Product number:	101.813-12 – including <i>Taq</i> polymerase 101.813-12u – without <i>Taq</i> polymerase
Lot number:	11K
Expiry date:	2012-May-01
Number of tests:	12
Number of wells per test:	2

CHANGES COMPARED TO THE PREVIOUS DRB1*14:01/14:54 Lot(72G).

The DRB1*14:01/14:54 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The DRB1*14:01 allele will give rise to a unique amplification pattern by the primers in the DRB1*14:01/14:54 kit; www.ebi.ac.uk/imgt/hla, 2010-April-01, release 3.0.0.

RESOLUTION IN DRB1*14:01/14:54 HOMOZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DRB1*14:01/14:54 SUBTYPINGS BY NON-DRB1*14:01/14:54 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.