

## **Olerup SSP<sup>®</sup> DRB1\*1401/1454**

<b>Product number:</b>	<b>101.813-12 – including <i>Taq</i> polymerase 101.813-12u – without <i>Taq</i> polymerase</b>
<b>Lot number:</b>	<b>72G</b>
<b>Expiry date:</b>	<b>2012-January-01</b>
<b>Number of tests:</b>	<b>12</b>
<b>Number of wells per test:</b>	<b>2</b>

**THE NUMBER OF WELLS** is 2 wells in an 8-well plate.

### **ALLELE COVERAGE:**

The DRB1\*1401 allele will give rise to a unique amplification pattern by the primers in the DRB1\*1401/1454 kit; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2009-October-19, release 2.27.0.

### **RESOLUTION IN DRB1\*1401/1454 HOMOZYGOTES:**

Excellent.

### **INFLUENCE ON THE INTERPRETATION OF DRB1\*1401/1454 SUBTYPINGS BY NON-DRB1\*1401/1454 ALLELES:**

None of importance.