

## **Olerup SSP<sup>®</sup> DRB1\*14:01/14:54**

**Product number:** 101.813-12 – including *Taq* polymerase  
101.813-12u – without *Taq* polymerase  
**Lot number:** 12L  
**Expiry date:** 2013-July-01  
**Number of tests:** 12  
**Number of wells per test:** 2

### **CHANGES COMPARED TO THE PREVIOUS DRB1\*14:01/14:54 Lot(11K).**

The DRB1\*14:01/14:54 primer set is unchanged compared to the previous lot.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

The DRB1\*14:01 allele will give rise to a unique amplification pattern by the primers in the DRB1\*14:01/14:54 kit; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2010-October-15, release 3.2.0.

### **RESOLUTION IN DRB1\*14:01/14:54 HOMOZYGOTES:**

Excellent.

### **INFLUENCE ON THE INTERPRETATION OF DRB1\*14:01/14:54 SUBTYPINGS BY NON-DRB1\*14:01/14:54 ALLELES:**

None of importance.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No suggestions received.