

Olerup SSP® DRB1*14:01/14:54

Product number:	101.813-12 – including <i>Taq</i> polymerase 101.813-12u – without <i>Taq</i> polymerase
Lot number:	97R
Expiry date:	2015-September-01
Number of tests:	12
Number of wells per test:	2

CHANGES COMPARED TO THE PREVIOUS DRB1*14:01/14:54 Lot(45N).

The DRB1*14:01/14:54 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The DRB1*14:01 and DRB1*14.54 alleles will give rise to unique amplification patterns by the primers in the DRB1*14:01/14:54 kit; www.ebi.ac.uk/imgt/hla, 2012-October-14, release 3.10.0.

RESOLUTION IN DRB1*14:01/14:54 HOMOZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DRB1*14:01/14:54 SUBTYPINGS BY NON-DRB1*14:01/14:54 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.