

101.813-12 – including *Taq* polymerase

101.813-12u – without *Taq* polymerase

Lot No.: **11X**

Olerup SSP® DRB1*14:01/14:54

Product number:	101.813-12 – including <i>Taq</i> polymerase 101.813-12u – without <i>Taq</i> polymerase
Lot number:	11X
Expiry date:	2017-April-01
Number of tests:	12
Number of wells per test:	2+1

CHANGES COMPARED TO THE PREVIOUS DRB1*14:01/14:54 Lot(97R).

Well	5'-primer	3'-primer	rationale
3	-	-	Negative Control.

THE NUMBER OF WELLS is increased from 2 to 3 wells.

ALLELE COVERAGE:

The DRB1*14:01 and DRB1*14.54 alleles will give rise to unique amplification patterns by the primers in the DRB1*14:01/14:54 kit;¹ www.ebi.ac.uk/imgt/hla, 2014-July-21, release 3.17.0.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN DRB1*14:01/14:54 HOMOZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DRB1*14:01/14:54 SUBTYPINGS BY NON-DRB1*14:01/14:54 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.