

101.813-12 – including *Taq* polymerase

101.813-12u – without *Taq* polymerase

Lot No.: **7D6**

**Olerup SSP® DRB1\*14:01/14:54**

**Product number:** 101.813-12 – including *Taq* polymerase  
 101.813-12u – without *Taq* polymerase  
**Lot number:** 7D6  
**Expiry date:** 2018-11-01  
**Number of tests:** 12  
**Number of wells per test:** 2+1

**CHANGES COMPARED TO THE PREVIOUS DRB1\*14:01/14:54 Lot(11X).**

Well	5'-primer	3'-primer	rationale
3	-	-	Updated negative control.

**THE NUMBER OF WELLS** is unchanged.

**ALLELE COVERAGE:**

The DRB1\*14:01 and DRB1\*14.54 alleles will give rise to unique amplification patterns by the primers in the DRB1\*14:01/14:54 kit;<sup>1</sup> [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2016-January-19, release 3.23.0.

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

**RESOLUTION IN DRB1\*14:01/14:54 HOMOZYGOTES:**

Excellent.

**INFLUENCE ON THE INTERPRETATION OF DRB1\*14:01/14:54 SUBTYPINGS BY NON-DRB1\*14:01/14:54 ALLELES:**

None of importance.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No suggestions received.