

101.813-12 – including *Taq* polymerase

101.813-12u – without *Taq* polymerase

Lot No.: **6H7**

Olerup SSP® DRB1*14:01/14:54

Product number:	101.813-12 – including <i>Taq</i> polymerase 101.813-12u – without <i>Taq</i> polymerase
Lot number:	6H7
Expiry date:	2023-04-01
Number of tests:	12
Number of wells per test:	2+1

CHANGES COMPARED TO THE PREVIOUS DRB1*14:01/14:54 LOT (0G7):

The DRB1*14:01/14:54 primer set is unchanged compared to the previous lot (**Lot No. 0G7**).

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The DRB1*14:01 and DRB1*14:54 alleles will give rise to unique amplification patterns by the primers in the DRB1*14:01/14:54 kit;¹ www.ebi.ac.uk/imgt/hla, 2019-January-23, release 3.35.0.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN DRB1*14:01/14:54 HOMOZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DRB1*14:01/14:54 SUBTYPINGS BY NON-DRB1*14:01/14:54 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.

Changes in revision R01 compared to R00:

1. The expiration date has been altered due to shelf-life extension.