

Release Note - Olerup SSP™ DRB3 – Lot No. Y16

Product number: 101.121-24/04 – licensed for PCR
101.121-24u/04u – not licensed for PCR
Lot number: Y16
Expiry date: 2009-July-01
Number of tests: 24 tests – Product No. 101.121-24
4 tests – Product No. 101.121-04
Number of tubes per test: 24

CHANGES COMPARED TO THE PREVIOUS DRB3 LOT (V61):

Tube	5'-primer	3'-primer	rationale
21	-	-	Changed to the 430 bp internal positive control primer pair to decrease primer dimer formation tendency.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DRB3*0101 to 0111, DRB3*0201 to 0222 and DRB3*0301 to 0303 alleles, i.e. all the currently recognized DRB3 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2007-July-09, release 2.18.0.

RESOLUTION IN DRB3 HOMOZYGOTES:

Very good.

INFLUENCE ON THE INTERPRETATION OF DRB3 SUBTYPINGS BY NON-DRB3 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.