

DRB4
101.122-24/06 – licensed for PCR
101.122-24u/06u – not licensed for PCR
Lot No.: **Y01**

1

Olerup SSP™ DRB4

Product number: 101.122-24/06 – licensed for PCR
101.122-24u/06u – not licensed for PCR
Lot number: Y01
Expiry date: 2009-June-01
Number of tests: 24 test – Product No. 101.122-24
6 tests – Product No. 101.122-06
Number of tubes per test: 12

CHANGES COMPARED TO THE PREVIOUS DRB4 LOT (V63):

Tube	5'-primer	3'-primer	rationale
1	-	Modified	Increased yield of specific PCR product.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DRB4*01010101 to DRB4*0301N, i.e. all the currently recognized DRB4 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2007-July-09, release 2.18.0.

RESOLUTION IN DRB4 HOMOZYGOTES:

The DRB4*01010101,01010101 and *01010101,0301N genotypes give rise to identical amplification patterns.

The DRB4*01010101,0103 and *0103,0301N genotypes give rise to identical amplification patterns.

INFLUENCE ON THE INTERPRETATION OF DRB4 SUBTYPINGS BY NON-DRB4 ALLELES:

The DRB1*150101 to DRB1*1522 and DRB1*160101 to DRB1*160502 and DRB1*1607 to DRB1*1611 alleles might be faintly amplified by primer mix 7.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.