

Olerup SSP[®] DRB4*01030102N

Product number:	101.811-12 – including <i>Taq</i> polymerase 101.811-12u – without <i>Taq</i> polymerase
Lot number:	41G
Expiry date:	2011-October-01
Number of tests:	12
Number of wells per test:	2

CHANGES COMPARED TO THE PREVIOUS DRB4*01030102N Lot (73E):

The DRB4*01030102N primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged. The plate layout has been changed to 2 wells in an 8-well plate.

ALLELE COVERAGE:

The DRB4*01030102N allele will give rise to a unique amplification pattern by the primers in the DRB4*01030102N kit; www.ebi.ac.uk/imgt/hla, 2009-July-17, release 2.26.0.

RESOLUTION IN DRB4*01030102N HOMOZYGOTES:

Excellent

INFLUENCE ON THE INTERPRETATION OF DRB4*01030102N SUBTYPINGS BY NON-DRB4*01030102N ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.