

101.811-12 – including *Taq* polymerase

101.811-12u – without *Taq* polymerase

Lot No.: 09V

Olerup SSP® DRB4*01:03:01:02N

Product number: 101.811-12 – including *Taq* polymerase
101.811-12u – without *Taq* polymerase

Lot number: 09V

Expiry date: 2016-May-01

Number of tests: 12

Number of wells per test: 2+1

CHANGES COMPARED TO THE PREVIOUS DRB4*01:03:01:02N LOT (09V):

Well	5'-primer	3'-primer	rationale
3	New	New	Negative Control.

THE NUMBER OF WELLS is increased from 2 to 3.

ALLELE COVERAGE:

The DRB4*01:03:01:02N allele will give rise to a unique amplification pattern by the primers in the DRB4*01:03:01:02N kit; www.ebi.ac.uk/imgt/hla, 2013-October-25, release 3.14.0.

RESOLUTION IN DRB4*01:03:01:02N HOMOZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DRB4*01:03:01:02N SUBTYPINGS BY NON-DRB4*01:03:01:02N ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.