

101.811-12 – including *Taq* polymerase101.811-12u – without *Taq* polymerase

Lot No.: 23Y

**Olerup SSP® DRB4\*01:03:01:02N**

**Product number:** 101.811-12 – including *Taq* polymerase  
101.811-12u – without *Taq* polymerase  
**Lot number:** 23Y  
**Expiry date:** 2017-September-01  
**Number of tests:** 12  
**Number of wells per test:** 2+1

**CHANGES COMPARED TO THE PREVIOUS DRB4\*01:03:01:02N LOT (09V):**

Well	5'-primer	3'-primer	rationale
3	-	-	Updated negative control.

**THE NUMBER OF WELLS** is unchanged.

**ALLELE COVERAGE:**

The DRB4\*01:03:01:02N allele will give rise to a unique amplification pattern by the primers in the DRB4\*01:03:01:02N kit; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2015-January-19, release 3.19.0<sup>1</sup>.

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

**RESOLUTION IN DRB4\*01:03:01:02N HOMOZYGOTES:**

Excellent.

**INFLUENCE ON THE INTERPRETATION OF DRB4\*01:03:01:02N SUBTYPINGS BY NON-DRB4\*01:03:01:02N ALLELES:**

None of importance.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.