

101.811-12 – including *Taq* polymerase

101.811-12u – without *Taq* polymerase

Lot No.: **6E4**

## **Olerup SSP® DRB4\*01:03:01:02N**

<b>Product number:</b>	101.811-12 – including <i>Taq</i> polymerase 101.811-12u – without <i>Taq</i> polymerase
<b>Lot number:</b>	6E4
<b>Expiry date:</b>	2019-06-01
<b>Number of tests:</b>	12
<b>Number of wells per test:</b>	2+1

The DRB4\*01:03:01:02N primer set is unchanged compared to the previous *Olerup SSP®* DRB4\*01:03:01:02N (Lot No. 23Y).

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

The DRB4\*01:03:01:02N allele will give rise to a unique amplification pattern by the primers in the DRB4\*01:03:01:02N kit; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2016-July-14, release 3.25.0<sup>1</sup>.

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

### **RESOLUTION IN DRB4\*01:03:01:02N HOMOZYGOTES:**

Excellent.

### **INFLUENCE ON THE INTERPRETATION OF DRB4\*01:03:01:02N SUBTYPINGS BY NON-DRB4\*01:03:01:02N ALLELES:**

None of importance.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.