

## **Olerup SSP<sup>®</sup> DRB5**

Product number:	101.123-24/06 – including <i>Taq</i> pol. 101.123-24u/06u – without <i>Taq</i> pol.
Lot number:	14M
Expiry date:	2013-October-01
Number of tests:	24 test – Product No. 101.123-24 6 tests – Product No. 101.123-06
Number of wells per test:	16

### **CHANGES COMPARED TO THE PREVIOUS DRB5 LOT (88G):**

The DRB5 primer set is unchanged compared to the previous lot.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

DRB5\*01:01 to DRB5\*01:14 and DRB5\*02:02 to DRB5\*02:05, i.e. all the currently recognized DRB5 alleles, give rise to unique amplification patterns, [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2011-January-14, release 3.3.0.

### **RESOLUTION IN DRB5 HOMOZYGOTES:**

Excellent.

### **INFLUENCE ON THE INTERPRETATION OF DRB5 SUBTYPINGS BY NON-DRB5 ALLELES:**

None frequently occurring.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.