

Olerup SSP[®] DRB5

| | |
|---------------------------|--|
| Product number: | 101.123-24/06 – including <i>Taq</i> pol. 101.123-24u/06u – without <i>Taq</i> pol. |
| Lot number: | 35N |
| Expiry date: | 2014-August-01 |
| Number of tests: | 24 test – Product No. 101.123-24/24u 6 tests – Product No. 101.123-06/06u |
| Number of wells per test: | 16 |

CHANGES COMPARED TO THE PREVIOUS DRB5 LOT (14M):

The DRB5 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DRB5*01:01 to DRB5*01:14 and DRB5*02:02 to DRB5*02:05, i.e. all the currently recognized DRB5 alleles, give rise to unique amplification patterns, www.ebi.ac.uk/imgt/hla, 2012-January-12, release 3.7.0.

RESOLUTION IN DRB5 HOMOZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DRB5 SUBTYPINGS BY NON-DRB5 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.