

101.812-12 – including *Taq* polymerase101.812-12u – without *Taq* polymeraseLot No.: **94K****Olerup SSP<sup>®</sup> DRB5\*01:08N**

Product number:	101.812-12 – including <i>Taq</i> polymerase 101.812-12u – without <i>Taq</i> polymerase
Lot number:	94K
Expiry date:	2013-July-01
Number of tests:	12
Number of wells per test:	2

**CHANGES COMPARED TO THE PREVIOUS DRB5\*01:08N Lot (42G):**

The DRB5\*01:08N primer set is unchanged compared to the previous lot.

**THE NUMBER OF WELLS** is unchanged.

**ALLELE COVERAGE:**

The DRB5\*01:08N allele will give rise to a unique amplification pattern by the primers in the DRB5\*01:08N kit; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2010-October-15, release 3.2.0.

**RESOLUTION IN DRB5\*01:08N HOMOZYGOTES:**

Excellent

**INFLUENCE ON THE INTERPRETATION OF DRB5\*01:08N SUBTYPINGS BY NON-DRB5\*01:08N ALLELES:**

None of importance.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.