

Olerup SSP[®] DRB5*01:08N

Product number:	101.812-12 – including <i>Taq</i> polymerase 101.812-12u – without <i>Taq</i> polymerase
Lot number:	44N
Expiry date:	2014-August-01
Number of tests:	12
Number of wells per test:	2

CHANGES COMPARED TO THE PREVIOUS DRB5*01:08N Lot (94K):

The DRB5*01:08N primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The DRB5*01:08N allele will give rise to a unique amplification pattern by the primers in the DRB5*01:08N kit; www.ebi.ac.uk/imgt/hla, 2012-January-12, release 3.7.0.

RESOLUTION IN DRB5*01:08N HOMOZYGOTES:

Excellent

INFLUENCE ON THE INTERPRETATION OF DRB5*01:08N SUBTYPINGS BY NON-DRB5*01:08N ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.