

101.812-12 – including *Taq* polymerase

101.812-12u – without *Taq* polymerase

Lot No.: **18Y**

Olerup SSP® DRB5*01:08N

Product number: 101.812-12 – including *Taq* polymerase
101.812-12u – without *Taq* polymerase
Lot number: 18Y
Expiry date: 2017-September-01
Number of tests: 12
Number of wells per test: 2+1

CHANGES COMPARED TO THE PREVIOUS DRB5*01:08N LOT (97V):

| Well | 5'-primer | 3'-primer | rationale |
|------|-----------|-----------|--------------------------|
| 3 | - | - | Updated negative control |

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The DRB5*01:08N allele will give rise to a unique amplification pattern by the primers in the DRB5*01:08N kit; www.ebi.ac.uk/imgt/hla, 2015-January-19, release 3.19.0.

RESOLUTION IN DRB5*01:08N HOMOZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DRB5*01:08N SUBTYPINGS BY NON-DRB5*01:08N ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.