

101.845-12 – including *Taq* polymerase101.845-12u – without *Taq* polymerase

Lot No.: 31Y

**Olerup SSP® HLA-A\*01 Add-on**

**Product number:** 101.845-12 – including *Taq* polymerase  
101.845-12u – without *Taq* polymerase

**Lot number:** 31Y

**Expiry date:** 2017-October-01

**Number of tests:** 12

**Number of wells per test:** 2+1

THE NUMBER OF WELLS is 3 in an 8 well plate.

**ALLELE COVERAGE:**

The HLA-A\*01:01:01:01G and A\*01:04N alleles give different patterns in the HLA-A\*01 Add-on kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2015-January-19, release 3.19.0.

The HLA-A\*01 add-on kit cannot distinguish the silent mutations in the A\*01:01:01:01-01:01:68 alleles.

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

**RESOLUTION IN HOMO- AND HETEROZYGOTES:**

Excellent.

**INFLUENCE ON THE INTERPRETATION OF HLA-A\*01 ADD-ON SUBTYPINGS BY NON-HLA-A\*01 ALLELES:**

None of importance.