

101.847-12 – including *Taq* polymerase101.847-12u – without *Taq* polymerase

Lot No.: 33Y

Olerup SSP® HLA-A*03 Add-on

Product number:	101.847-12 – including <i>Taq</i> polymerase 101.847-12u – without <i>Taq</i> polymerase
Lot number:	33Y
Expiry date:	2017-November-01
Number of tests:	12
Number of wells per test:	2+1

THE NUMBER OF WELLS is 3 in an 8 well plate.

ALLELE COVERAGE:

The HLA-A*03:01 and A*03:21N alleles give different patterns in the HLA-A*03 Add-on kit¹; www.ebi.ac.uk/imgt/hla, 2015-January-19, release 3.19.0.

The HLA-A*03 add-on kit cannot distinguish the silent mutations in the A*03:01:01:01-03:01:29, 03:01:31-03:01:47, 03:01:51-03:01:53 alleles.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-A*03 ADD-ON SUBTYPINGS BY NON-HLA-A*03 ALLELES:

None of importance.