

101.411-24/06 – including *Taq* polymerase101.411-24/06u – without *Taq* polymeraseLot No.: **94E****Olerup SSP[®] HLA-A*01**

Product number: 101.411-24/06– including *Taq* pol.
101.411-24u/06u – without *Taq* pol.

Lot number: 94E

Expiry date: 2010-August-01

Number of tests: 24 test – Product No. 101.411-24
6 tests – Product No. 101.411-06

Number of wells per test: 24

CHANGES COMPARED TO THE PREVIOUS HLA-A*01 LOT (X98):

Well	5'-primer	3'-primer	rationale
8	-	Added	Primer added for the A*01010102N allele.
10	Added	Added	Primer pair added for the A*0129 allele.
11	Added	Added	Primer pair added for the A*0128 allele.
19	Added	Added	Primer pair added for the A*0130 allele.
23	Added	Added	Primer pair added for the A*0130 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

A*0101 to A*0130, i.e. all the currently recognized HLA-A*01 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-July-11, release 2.22.0.

RESOLUTION IN HLA-A*01 HOMO- AND HETEROZYGOTES:

The A*0101,0101, A*0101,0102 and A*0102,0102 genotypes give rise to unique amplification patterns.

INFLUENCE ON THE INTERPRETATION OF HLA-A*01 SUBTYPINGS BY NON-HLA-A*01 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.