

Lot No.: **30F**

Olerup SSP[®] HLA-A*11

Product number: 101.416-12/04 – including *Taq* pol.
101.416-12u/04u – without *Taq* pol.
Lot number: 30F
Expiry date: 2010-November-01
Number of tests: 12 tests – Product No. 101.416-12/12u
4 tests – Product No. 101.416-04/04u
Number of Wells per test: 32

CHANGES COMPARED TO THE PREVIOUS HLA-A*11 LOT (Y15):

Well	5'-primer	3'-primer	rationale
3	-	Added	Primer added for the A*1133 allele.
8	Added	-	Primer added for the A*1134 allele.
13	Added	Added	Primer pair added for the A*1134 allele.
15	Added	Added	Primer pair added for the A*112402 allele.
13	Added	Added	Primer pair added for the A*1136 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

A*1101 to A*1136, i.e. all the currently recognized HLA-A*11 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-July-11, release 2.22.0.

RESOLUTION IN HLA-A*11 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-A*11 SUBTYPINGS BY NON-HLA-A*11 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.