

## **Olerup SSP<sup>®</sup> HLA-A\*11 Add-on**

<b>Product number:</b>	<b>101.842-12 – including <i>Taq</i> polymerase</b> <b>101.842-12u – without <i>Taq</i> polymerase</b>
<b>Lot number:</b>	<b>21N</b>
<b>Expiry date:</b>	<b>2014-August-01</b>
<b>Number of tests:</b>	<b>12</b>
<b>Number of wells per test:</b>	<b>3</b>

**THE NUMBER OF WELLS** is 3 in an 8 well plate.

### **ALLELE COVERAGE:**

The HLA-A\*11:02, 11:77 and 11:110 alleles give different patterns in the HLA-A\*11 Add-on kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2012-January-12, release 3.7.0.

<sup>1</sup>The HLA-A\*11 add-on kit cannot distinguish the silent mutations in the A\*11:02:01-11:02:03 alleles.

### **RESOLUTION IN HOMO- AND HETEROZYGOTES:**

Excellent.

### **INFLUENCE ON THE INTERPRETATION OF HLA-A\*11 ADD-ON SUBTYPINGS BY NON-HLA-A\*11 ALLELES:**

None of importance.