

Olerup SSP[®] HLA-A*23 Add-on

Product number:	101.843-12 – including <i>Taq</i> polymerase 101.843-12u – without <i>Taq</i> polymerase
Lot number:	55N
Expiry date:	2014-August-01
Number of tests:	12
Number of wells per test:	3

THE NUMBER OF WELLS is 3 in an 8 well plate.

ALLELE COVERAGE:

The HLA-A*23:01 and 23:17 alleles give different patterns in the HLA-A*23 Add-on kit¹; www.ebi.ac.uk/imgt/hla, 2012-January-12, release 3.7.0.

¹The HLA-A*23 add-on kit cannot distinguish the silent mutations in the A*23:01:01-23:01:07 alleles.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-A*23 ADD-ON SUBTYPINGS BY NON-HLA-A*23 ALLELES:

None of importance.