

101.843-12 – including *Taq* polymerase

101.843-12u – without *Taq* polymerase

Lot No.: **70X**

Olerup SSP® HLA-A*23 Add-on

Product number:	101.843-12 – including <i>Taq</i> polymerase 101.843-12u – without <i>Taq</i> polymerase
Lot number:	70X
Expiry date:	2017-May-01
Number of tests:	12
Number of wells per test:	3+1

CHANGES COMPARED TO THE PREVIOUS HLA-A*23 ADD-ON LOT (48S):

Well	5'-primer	3'-primer	rationale
4	-	-	Negative Control.

THE NUMBER OF WELLS is increased from 3 to 4 wells.

ALLELE COVERAGE:

The HLA-A*23:01 and 23:17 alleles give different patterns in the HLA-A*23 Add-on kit¹; www.ebi.ac.uk/imgt/hla, 2014-July-25, release 3.17.0.

¹The HLA-A*23 add-on kit cannot distinguish the silent mutations in the A*23:01:01-23:01:14 alleles.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-A*23 ADD-ON SUBTYPINGS BY NON-HLA-A*23 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.