

101.843-12 – including *Taq* polymerase101.843-12u – without *Taq* polymeraseLot No.: **6E6****Olerup SSP® HLA-A\*23 Add-on**

**Product number:** 101.843-12 – including *Taq* polymerase  
101.843-12u – without *Taq* polymerase

**Lot number:** 6E6

**Expiry date:** 2019-06-01

**Number of tests:** 12

**Number of wells per test:** 3+1

The HLA-A\*23 Add-on primer set is unchanged compared to the previous *Olerup SSP® HLA-A\*23 Add-on (Lot No. 8D6)*.

**THE NUMBER OF WELLS** is unchanged.

**ALLELE COVERAGE:**

The HLA-A\*23:01 and 23:17 alleles give different patterns in the HLA-A\*23 Add-on kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2016-July-14, release 3.25.0.

<sup>1</sup>The HLA-A\*23 add-on kit cannot distinguish the silent mutations in the A\*23:01:01-23:01:15 and 23:01:17-23:01:19 alleles.

**RESOLUTION IN HOMO- AND HETEROZYGOTES:**

Excellent.

**INFLUENCE ON THE INTERPRETATION OF HLA-A\*23 ADD-ON SUBTYPINGS BY NON-HLA-A\*23 ALLELES:**

None of importance.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.