

101.843-12 – including *Taq* polymerase

101.843-12u – without *Taq* polymerase

Lot No.: **8D6**

## **Olerup SSP® HLA-A\*23 Add-on**

**Product number:** 101.843-12 – including *Taq* polymerase  
101.843-12u – without *Taq* polymerase

**Lot number:** 8D6

**Expiry date:** 2018-09-01

**Number of tests:** 12

**Number of wells per test:** 3+1

### **CHANGES COMPARED TO THE PREVIOUS HLA-A\*23 ADD-ON LOT (70X):**

Well	5'-primer	3'-primer	rationale
4	-	-	Updated negative control.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

The HLA-A\*23:01 and 23:17 alleles give different patterns in the HLA-A\*23 Add-on kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2016-January-19, release 3.23.0.

<sup>1</sup>The HLA-A\*23 add-on kit cannot distinguish the silent mutations in the A\*23:01:01-23:01:15 and 23:01:17-23:01:18 alleles.

### **RESOLUTION IN HOMO- AND HETEROZYGOTES:**

Excellent.

### **INFLUENCE ON THE INTERPRETATION OF HLA-A\*23 ADD-ON SUBTYPINGS BY NON-HLA-A\*23 ALLELES:**

None of importance.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.