

101.846-12 – including *Taq* polymerase101.846-12u – without *Taq* polymerase

Lot No.: 1E4

**Olerup SSP® HLA-A\*24 Add-on**

**Product number:** 101.846-12 – including *Taq* polymerase  
101.846-12u – without *Taq* polymerase

**Lot number:** 1E4

**Expiry date:** 2018-12-01

**Number of tests:** 12

**Number of wells per test:** 3+1

This lot is manufactured using white plastic trays.

**THE NUMBER OF WELLS** is 4 in an 8 well plate.

**ALLELE COVERAGE:**

The HLA-A\*24:09N, A\*24:11N and A\*24:02:01G alleles give different patterns in the HLA-A\*24 Add-on kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2016-January-19, release 3.23.0.

The HLA-A\*24 add-on kit cannot distinguish the silent mutations in the A\*24:02:01:01-24:02:96 alleles.

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

**RESOLUTION IN HOMO- AND HETEROZYGOTES:**

Excellent.

**INFLUENCE ON THE INTERPRETATION OF HLA-A\*24 ADD-ON SUBTYPINGS BY NON-HLA-A\*24 ALLELES:**

None of importance.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

In primer mix 2, the primer pair was exchanged for the A\*24:11N allele.