

Lot No.: **8E0****Olerup SSP® HLA-A\*24 Add-on**

**Product number:** 101.846-12 – including *Taq* polymerase  
101.846-12u – without *Taq* polymerase

**Lot number:** 8E0

**Expiry date:** 2019-05-01

**Number of tests:** 12

**Number of wells per test:** 3+1

**CHANGES COMPARED TO THE PREVIOUS HLA-A\*24 ADD-ON LOT (1E4):**

Well	5'-primer	3'-primer	rationale
1	Modified	-	5'-primer modified for improved HLA-specific amplification.

**THE NUMBER OF WELLS** is unchanged.

**ALLELE COVERAGE:**

The HLA-A\*24:09N, A\*24:11N and A\*24:02:01G alleles give different patterns in the HLA-A\*24 Add-on kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2016-July-14, release 3.25.0.

The HLA-A\*24 add-on kit cannot distinguish the silent mutations in the A\*24:02:01:01-24:02:97 alleles.

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

**RESOLUTION IN HOMO- AND HETEROZYGOTES:**

Excellent.

**INFLUENCE ON THE INTERPRETATION OF HLA-A\*24 ADD-ON SUBTYPINGS BY NON-HLA-A\*24 ALLELES:**

None of importance.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.