

101.846-12 – including *Taq* polymerase

101.846-12u – without *Taq* polymerase

Lot No.: 3F2

Olerup SSP[®] HLA-A*24 Add-on

Product number:	101.846-12 – including <i>Taq</i> polymerase 101.846-12u – without <i>Taq</i> polymerase
Lot number:	3F2
Expiry date:	2021-07-01
Number of tests:	12
Number of wells per test:	3+1

The HLA-A*24 Add-on primer set is unchanged compared to the previous *Olerup SSP[®] HLA-A*24 Add-on (Lot No. 8E0)*.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The HLA-A*24:09N, A*24:11N and A*24:02:01G alleles give different patterns in the HLA-A*24 Add-on kit¹; www.ebi.ac.uk/imgt/hla, 2016-October-14, release 3.26.0.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-A*24 ADD-ON SUBTYPINGS BY NON-HLA-A*24 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.

Changes in revision R01 compared to R00:

1. The expiration date has been altered due to extension of shelf-life.