

Olerup SSP[®] HLA-A*2409N

Product number: 101.841-12 – including *Taq* polymerase
101.841-12u – without *Taq* polymerase
Lot number: 27G
Expiry date: 2011-August-01
Number of tests: 12
Number of wells per test: 2

CHANGES COMPARED TO THE PREVIOUS HLA-A*2409N LOT (79E):
The HLA-A*2409N primer set is unchanged.

The plate layout of the HLA-A*2409N kit has been changed to 2 wells in an 8-well plate.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The HLA-A*2409N allele will give rise to a unique amplification pattern by the primers in the HLA-A*2409N kit; www.ebi.ac.uk/imgt/hla, 2009-July-17, release 2.26.0.

RESOLUTION IN HLA-A*2409N HOMO- AND HETEROZYGOTES:
Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-A*2409N SUBTYPINGS BY NON-HLA-A*2409N ALLELES:
None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:
No suggestions received.