

Olerup SSP[®] HLA-A*24:09N

Product number:	101.841-12 – including <i>Taq</i> polymerase 101.841-12u – without <i>Taq</i> polymerase
Lot number:	98R
Expiry date:	2015-September-01
Number of tests:	12
Number of wells per test:	2

CHANGES COMPARED TO THE PREVIOUS HLA-A*24:09N LOT (28R):
The HLA-A*24:09N primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The HLA-A*24:09N allele will give rise to a unique amplification pattern by the primers in the HLA-A*24:09N kit; www.ebi.ac.uk/imgt/hla, 2012-October-14, release 3.10.0.

RESOLUTION IN HLA-A*24:09N HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-A*24:09N SUBTYPINGS BY NON-HLA-A*24:09N ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.