

101.841-12 – including *Taq* polymerase

101.841-12u – without *Taq* polymerase

Lot No.: **48V**

Olerup SSP® HLA-A*24:09N

Product number:	101.841-12 – including <i>Taq</i> polymerase 101.841-12u – without <i>Taq</i> polymerase
Lot number:	48V
Expiry date:	2016-August-01
Number of tests:	12
Number of wells per test:	2+1

CHANGES COMPARED TO THE PREVIOUS HLA-A*24:09N LOT (98R):

Well	5'-primer	3'-primer	rationale
3	New	New	Negative Control.

THE NUMBER OF WELLS is increased from 2 to 3 wells.

ALLELE COVERAGE:

The HLA-A*24:09N allele will give rise to a unique amplification pattern by the primers in the HLA-A*24:09N kit; www.ebi.ac.uk/imgt/hla, 2014-January-17¹, release 3.15.0.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HLA-A*24:09N HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-A*24:09N SUBTYPINGS BY NON-HLA-A*24:09N ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.