

HLA-A*24:09N
101.841-12 – including *Taq* polymerase
101.841-12u – without *Taq* polymerase
Lot No.: 0H7

Olerup SSP[®] HLA-A*24:09N

Product number:	101.841-12 – including <i>Taq</i> polymerase 101.841-12u – without <i>Taq</i> polymerase
Lot number:	0H7
Expiry date:	2021-04-01
Number of tests:	12
Number of wells per test:	2+1

The HLA-A*24:09N primer set is unchanged compared to the previous *Olerup SSP[®]* HLA-A*24:09N (Lot No. 2E4).

THE NUMBER OF WELLS is 3 in an 8 well plate.

ALLELE COVERAGE:

The HLA-A*24:09N allele will give rise to a unique amplification pattern by the primers in the HLA-A*24:09N kit; www.ebi.ac.uk/imgt/hla, 2018-July-11¹, release 3.33.0.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HLA-A*24:09N HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-A*24:09N SUBTYPINGS BY NON-HLA-A*24:09N ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.