

Lot No.: **8H6**

Olerup SSP[®] HLA-A*24:09N

Product number: 101.841-12 – including *Taq* polymerase
101.841-12u – without *Taq* polymerase
Lot number: 8H6
Expiry date: 2023-06-01
Number of tests: 12
Number of wells per test: 3+1

CHANGES COMPARED TO THE PREVIOUS HLA-A*24:09N LOT (0H7):

Well	5'-primer	3'-primer	rationale
3	-	-	Negative control moved to well 4. Primer pair added for improved resolution of A*24 alleles.
4	-	-	Negative control added from well 3.

THE NUMBER OF WELLS is 4 in an 8 well plate.

ALLELE COVERAGE:

The HLA-A*24:09N allele will give rise to a unique amplification pattern by the primers in the HLA-A*24:09N kit; www.ebi.ac.uk/imgt/hla, 2019-January-23¹, release 3.35.0.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HLA-A*24:09N HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-A*24:09N SUBTYPINGS BY NON-HLA-A*24:09N ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.

Changes in revision R01 compared to R00:

1. The expiration date has been altered due to shelf-life extension.