

Olerup SSP[®] HLA-A*25

Product number: 101.423-06 – including *Taq* polymerase
101.423-06u – without *Taq* polymerase
Lot number: 23F
Expiry date: 2010-November-01
Number of tests: 6
Number of Wells per test: 8

CHANGES COMPARED TO THE PREVIOUS HLA-A*25 LOT (15E):

Well	5'-primer	3'-primer	rationale
7	Modified	Modified	Increased specificity of specific primer pair.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

A*2501 to A*2506, i.e. all the currently recognized HLA-A*25 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-October-10, release 2.23.0.

RESOLUTION IN HLA-A*25 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-A*25 SUBTYPINGS BY NON-HLA-A*25 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.