

Lot No.: **20G**

Olerup SSP[®] HLA-A*25

Product number: 101.423-06 – including *Taq* polymerase
101.423-06u – without *Taq* polymerase
Lot number: 20G
Expiry date: 2011-September-01
Number of tests: 6
Number of wells per test: 10

CHANGES COMPARED TO THE PREVIOUS HLA-A*25 LOT (23F):

Well	5'-primer	3'-primer	rationale
9	New	New	New primer pair for the A*2508 allele.
10	New	New	New primer pair for the A*2507 allele.

THE NUMBER OF WELLS is increased from 8 to 10, and the plate layout is changed from an 8 well plate to a 16 well plate.

ALLELE COVERAGE:

A*2501 to A*2508, i.e. all the currently recognized HLA-A*25 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2009-July-17, release 2.26.0.

RESOLUTION IN HLA-A*25 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-A*25 SUBTYPINGS BY NON-HLA-A*25 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.