

## **Olerup SSP<sup>®</sup> HLA-A\*25**

**Product number:** 101.423-06 – including *Taq* polymerase  
101.423-06u – without *Taq* polymerase  
**Lot number:** 15E  
**Expiry date:** 2010-February-01  
**Number of tests:** 6  
**Number of Wells per test:** 8

### **CHANGES COMPARED TO THE PREVIOUS HLA-A\*25 LOT (X36):**

The HLA-A\*25 primer set is unchanged.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

A\*2501 to A\*2506, i.e. all the currently recognized HLA-A\*25 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2008-January-11, release 2.20.0.

### **RESOLUTION IN HLA-A\*25 HOMO- AND HETEROZYGOTES:**

Good.

### **INFLUENCE ON THE INTERPRETATION OF HLA-A\*25 SUBTYPINGS BY NON-HLA-A\*25 ALLELES:**

None of importance.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No modifications made.