

HLA-A\*29  
101.428-12 – licensed for PCR  
101.428-12u – not licensed for PCR  
Lot No.: **X81**

1

## Release Note - *Olerup* SSP™ HLA-A\*29 – Lot No. **X81**

Product number: 101.428-12 – licensed for PCR  
101.428-12u – not licensed for PCR  
Lot number: X81  
Expiry date: 2008-April-01  
Number of tests: 12  
Number of tubes per test: 16

### CHANGES COMPARED TO THE PREVIOUS HLA-A\*29 LOT (V56):

Tube	5'-primer	3'-primer	rationale
9	Added	Added	Primer pair added for the A*2916 allele.
16	New	New	New primer pair for A*2915 allele.

**THE NUMBER OF WELLS** has been increased from 15 to 16.

### ALLELE COVERAGE:

A\*2901 to A\*2916, i.e. all the currently recognized HLA-A\*29 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2007-April-12, release 2.17.0.

### RESOLUTION IN HLA-A\*29 HOMO- AND HETEROZYGOTES:

Excellent.

### INFLUENCE ON THE INTERPRETATION OF HLA-A\*29 SUBTYPINGS BY NON-HLA-A\*29 ALLELES:

None frequently occurring.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.