

Olerup SSP[®] HLA-A*29

Product number: 101.428-12 – including *Taq* polymerase
101.428-12u – without *Taq* polymerase
Lot number: 64E
Expiry date: 2010-May-01
Number of tests: 12
Number of wells per test: 16

CHANGES COMPARED TO THE PREVIOUS HLA-A*29 LOT (X81):

Well	5'-primer	3'-primer	rationale
8	Modified	-	Increased yield of specific PCR product.
10	Added	Added	Primer pair added for the A*2917 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

A*2901 to A*2917, i.e. all the currently recognized HLA-A*29 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-April-08, release 2.21.0.

RESOLUTION IN HLA-A*29 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-A*29 SUBTYPINGS BY NON-HLA-A*29 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.